


Simplicity+

A Minimum swivel radius	655mm
– with arm spacers	678mm
– with extended arm spacers	695mm
B Minimum wall to seat back	65mm
C Seat depth	360mm
– with extended seat	425mm
D Minimum folded width to footplate	370mm
E Minimum open width to edge of footplate	570mm
– large footplate (perch seat)	630mm
F Armrest width – external	550mm
– with arm spacers	590mm
– with extended arm spacers	637mm
G Armrest width – internal	460mm
– with arm spacers	516mm
– with extended arm spacers	562mm
H Seat back height	453mm
I Footplate to seat height	
– standard seat option	430-540mm
– perch seat option	677-877mm
J Minimum footplate height	70mm
Bottom parking space required	575mm (50°) to 690 at shallow angle (30°)
Minimum track intrusion into staircase	185mm
Minimum track overhang at top of stairs (typical at 45°)	150mm
Minimum track intrusion at bottom of stairs	250mm
Maximum weight limit	140kg/22 stone
– Hinge/Slide Track	120kg/19 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	100-240v (50-60Hz)
Staircase angle	27-50°



Minivator Stairlifts



User support: 08700 118282

Handicare, 82 First Avenue, Pensnett Estate, Kingswinford, West Midlands DY6 7FJ, United Kingdom
Textphone: (180015) 020 7490 8913 ext 05 Voicephone: (18002) 020 7490 8913 ext 05
Global Dealer Support: +44 (0) 844 225 3121 enquiries@handicare.co.uk www.handicare.co.uk

Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Minivator Group reserves the right to alter designs and specifications without prior notice.



Simplicity+
Straight stairlift



The benefits of a Minivator Simplicity+ Stairlift from Handicare

If you are thinking about installing a stairlift the Simplicity+ offers you a safe and reliable solution.

Handicare's Minivator Stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Simplicity+ straight stairlift offers you all the basic functions that you would expect plus the option to upgrade to a powered seat swivel and/or folding footplate.

If you have an obstruction at the base of the stairs, you can choose from a manually operated hinged track or the unique Slide Track that is only available from Handicare.

With a choice of six upholstery colours the Simplicity+ will compliment most colour schemes.

The track is fitted to the stairs, not the wall, so it is unlikely you will need to redecorate following the installation of your Simplicity+ stairlift.

Like all Handicare stairlifts the Simplicity+ has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the Simplicity+ also comes with two remote controls which means, once folded, you can move it along the track freeing up space in your hallway.



Safety As Standard

Your Simplicity+ stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Simplicity+ stairlift has an easy release retractable seat belt for added security and peace of mind. If you prefer you can have an easy fasten Velcro seat belt.

Easy to Use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Simplicity+ will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.